

Dental Chair Test

1. What is the dental compressor air source pressure?
 - a. 30-50 psi
 - b. 40-50 psi
 - c. 60-80 psi
 - d. 60-100 psi
2. What requires replacement if it becomes clogged and pressure drop across the dental unit exceeds 10psi differential pressure. ?
 - a. Air filter
 - b. Water filter
 - c. Pressure transducer
 - d. Flush valve
3. What is used during the lubrication of dental handpiece components?
 - a. Thread lock
 - b. Vaseline
 - c. Lithium grease
 - d. Silicone
4. What instruments can be used to remove springs and o-rings from a 3-way syringe?
 - a. Dental pick
 - b. Paper clip
 - c. A and B
 - d. Screw driver
5. What component powers the hydraulic motor within the chair?
 - a. Solenoid
 - b. Data communication System
 - c. Capacitor
 - d. Potentiometers
6. What component senses and interprets when the chair is being raised or lowered within the chair?
 - a. Solenoid
 - b. Data communication System
 - c. Capacitor
 - d. Potentiometers
7. What component is considered “the brainbox” within the chair?

- a. Solenoid
- b. Data communication System
- c. Capacitor
- d. Potentiometers

8. SCENARIO: A dentist calls you on the phone. She/He tells you that the chair is not working and they're constantly reprogramming the memory on the touchpad. What is most likely NOT the cause of this problem?

- a. PCB is defective
- b. User error. Chair not correctly programmed
- c. Defective solenoid
- d. Defective base or back pot

9. SCENARIO: A dentist calls you on the phone. She/He tells you that it won't go upward. What is most likely NOT the cause of this problem?

- a. Hydrostatic lock condition
- b. User error. Chair not correctly programmed
- c. Defective solenoid
- d. Chair not correctly programmed

10. A dentist calls you on the phone. She/He tells you that it won't go up. What is most likely NOT the cause of this problem?

- a. Bad touch pad or footswitch
- b. Stop plate activated
- c. Chair not correctly programmed
- d. Defective motor capacitor

Answer Key

1. c
2. a
3. d
4. c
5. c
6. d
7. b
8. c
9. a
10. c